



InnoScreen

6 Panel Drug Screen

Verified to AS/NZ 4308:2008



Verified to Appendix B of the AS/NZ4308:2008 Drug Testing Standard, InnoScreen offers users the collector activated key (split cup) for convenience with its unique lid which the donor snaps shut before returning to the collector. Accurate and cost effective 6 panel drug screening- providing results in minutes, with no spills, no fuss, no mess.

Simply collect sample, activate with key and read temperature and integrity results before viewing test results

CONVENIENCE

- Tests for the six drugs of abuse as per AS/NZ4308:2008 requirements/cut-offs
- Incorporates collection and testing
- Temperature strip on cup
- Adulterant strip included to defeat specimen adulteration
- Simple procedure, collect sample read test in 5 minutes, no timing required
- Unique lid requires collector activation

PERFORMANCE

- Each test has procedural control built in Colloidal gold technology
- >99% correlation to GC/MS at 97.5% confidence level
- Tested for cross-reactivity against 500 drugs
- **AS/NZ Standard 4308:2008 cut-offs**
- **Verified under Appendix B to the AS/NZ 4308:2008 Standard**

CUT-OFF LEVELS

- Cut-off levels to AS/NZ Standard 4308:2008 for immunoassay testing : -
Amphetamines/Methamphetamines 300ng/ml,
Benzodiazepines 200ng/ml, Cocaine 300ng/ml
Opiates 300ng/ml, THC (Marijuana) 50ng/ml

AMP

amphetamines 300ng/ml

BZO

benzodiazepines 200ng/ml

COC

cocaine 300ng/ml

METH

methamphetamine 300ng/ml

OPI

opiates 300ng/ml

THC

marijuana 50ng/ml

Plus 3 Integrity Tests

- **Oxidants**
- **Creatinine**
- **pH (alkaline & acidity)**



Innoscreen

6 Panel Drug Screen Verified to AS/NZ 4308:2008



Innoscreen 6 panel drug screen is a rapid urine screening immunoassay device for the detection of the drugs of abuse. Easy to use, accurate and cost effective with results available in minutes ensuring no messy spills.



Easy to Use


Open foil pouch, remove key from lid of device and give specimen cup to donor

1. Donor provides specimen returning sample to collector sealed with lid
2. Collector checks temperature strip is in range (32-38C)
3. Collector then activates with key
4. Collector then peels label from cup.
5. View results, checking integrity tests with colour card provided
6. Interpret drug screening results:
 - Control line (C) should appear in all test windows.
 - A second line (T) should then appear for all test windows (**negative**)

 - A positive/**non-negative** is indicated by the **ABSENCE** of a test line below the corresponding drug
 - line intensity (colour) is **NOT** important)



Reading Results



C = Control
T = Test

All positive results are presumptive & should be **confirmed by lab GC/MS (Gas Chromatograph/Mass Spectrometry)**

Line Intensity is **NOT** important

- Negative results can be read as soon as they appear
- Read non-negative results at 5 minutes
- Line intensity is **NOT** important
- All non-negative results are presumptive and should be confirmed by GC/MS (Gas Chromatograph/Mass Spectrometry)

Negative = 2 lines in each window
NON Negative = only a Control line visible in window