



SureStep EZI Cup™



Verified to Appdx B AS/NZ 4308:2008 6 Panel Split Cup Drug Screen

With no fuss, no mess and a hinged lid clicking shut - ensuring no spills, EZI Cup™ is Verified to Appendix B of the AS/NZ4308:2008 Drug Testing Standard. Offering convenience, accuracy and cost effective drug screening at your fingertips, EZI Cup™ is simple to use with negative results available in minutes .

CONVENIENCE

- Tests for the six drugs of abuse as per AS/NZ4308:2008 requirements/cut-offs
- Incorporates collection and testing
- Temperature strip on cup
- Adulterant strip mandatory to defeat specimen adulteration
- Simple procedure, collect sample, click shut.
- Collector retains complete control, activating cup with key
- Read test in 5 minutes, no timing required

PERFORMANCE

- Each test has procedural control built in Colloidal gold technology
- Tested for cross-reactivity against 500 drugs
- AS/NZ Standard 4308:2008 cut-offs
- Verified under Appendix B to the AS/NZ 4308:2008 Standard
- >99% correlation to GC/MS at 97.5% confidence level

CUT-OFF LEVELS

- Cut-off levels to AS/NZ Standard 4308:2008 for immunoassay testing : -
Amphetamines/Methamphetamines 300ng/ml,
Benzodiazepines 200ng/ml, Cocaine 300ng/ml
Opiates 300ng/ml, THC (Marijuana) 50ng/ml

AMP

amphetamines 300ng/ml

BZO

benzodiazepines 200ng/ml

COC

cocaine 300ng/ml

METH

methamphetamine 300ng/ml

OPI

opiates 300ng/ml

THC

marijuana 50ng/ml

Plus 6 Integrity Tests

- Creatinine
- Nitrites
- pH (alkaline & acidity)
- Bleach- (Oxidants)
- SG (Specific Gravity)
- Glutaraldehyde

SureStep EZI Cup™

Verified to Appdx B
AS/NZ 4308:2008
6 Panel Split Cup
Drug Screen



Perfect for workplace drug screening, user friendly, accurate and cost effective with results available in minutes, the EZI Cup™ 6 panel drug screen is a rapid urine screening immunoassay device for the detection of the drugs of abuse in urine.

User Friendly

Open foil pouch, **remove key from lid and retain**, - give specimen cup to donor











1. Donor provides specimen and returns sample to collector
2. Collector then checks hinged lid securely closed
3. Collector checks temperature strip is in range (32-38C)
4. Insert key into side of test device and push to activate
5. Peel off label, wait 3-5 minutes for results to appear
6. View results, check integrity tests with colour card provided
7. Interpret drug screening results:
 - Control line (C) should appear in all test windows.
 - A second line (T) should then appear for all test windows (**negative**)



- A positive/**non-negative** is indicated by the **ABSENCE** of a test line (T) below the corresponding drug - line intensity (colour) is not important

Reading Results

- Negative results can be read as soon as they appear
- Read non-negative results at 5 minutes
- Line intensity is **NOT** important
- All non-negative results are presumptive and should be confirmed by GC/MS (Gas Chromatograph/Mass Spectrometry)

C			C		C = Control T = Test
T			T		
NEGATIVE		POSITIVE		All positive results are presumptive & should be confirmed by lab GC/MS (Gas Chromatograph/Mass Spectrometry)	
C			INVALID		
T					

Line Intensity is **NOT** important

- **Negative** = 2 lines in each window
- **Window 3 : Must show 4 lines** (including Control Line) Negative for all drugs
- **NON Negative** = only a Control line visible in window and absence of Test line/s