

## PROSCREEN USERS

**Currently** we are unable to supply this urine drug test cup, as is the case for all suppliers of this product worldwide. The production of Proscreen/Microscreen urine drug testing cups worldwide has been dramatically impacted by a recent buy-out, causing the continued supply of this most widely used drug testing device to effectively "grind to a halt."

This in-turn has impacted ours and others ability to supply. The issue of supply or (non-supply) we are told by the manufacturers, is likely to extend into February 2020.

---

In the interim we are able to offer alternative drug screening options to fill the temporary gap due to the supply chain being halted - until this issue is resolved.

There are three distinct urine drug testing options available, all verified to Appdx B AS/NZ 4308:2008 and either tilt cup as per Proscreen, EcoCup II (direct delivery of the specimen) or the split cup (EZI Cup and Innoscreen.) -

### **1. Innoscreen Rapid Cup verified to AS/NZ 4308:2008.** 6 panel urine c/w 6 integrity tests

Based on the same principle as the Proscreen Cup :  
Collector screws the lid on and test is activated by tilting the cup onto the legs of the lid. Used by NSW Police exclusively for their urine drug testing regime. This option would provide a seamless processing option

### **2. ECOcup II** urine drug test cup 6 panel. Again verified to AS/NZS4308:2008

This cup is a good alternative, however the test process isn't performed the same as a tilt cup type. The Ecocup will start to process the test (integrity and drug test) immediately the sample is delivered into the cup.  
- Results are read on the label around the outer of the cup.



### **3. InnoScreen Split Cup verified to AS/NZ 4308:2008**

Offers the same performance as the SureStep EZI Cup or Medix Prosplit cup but has the added advantage of each drug class having its own test window.



### **4. Surestep EZI Cup - verified too AS/NZ 4308:2008**

Split cup 6 panel urine screen with a push button key used to deliver sufficient sample into the testing chamber (split cup) thus activating the control for the screening process



We are hoping this supply issue will be shortlived and you will return to your Proscreen use when it is finally sorted. Until then, Drugchek apologize for any inconvenience caused and offer as always, their full ongoing customer support and service during this transition time and are always happy to discuss any questions or concerns you may have.